

Very small visual appearance and non medical look, low volume design that makes the therapy easier for patients. Mini J5 is for patients that require regular therapy in a smaller and power efficient platform.



5 LPM Oxygen

%93 Oxygen Purity

330 Watts (Economy)

15.1 Kilograms (Light)

- Low power consumption
- Less heat production
- Silent operation

Technical specifications

Part Number	M0000025
Operating Voltage	220-240 V/ AC
Main Frequency	50-60 Hz
Power	< 310 W
Oxygen Flow	0.5-5 LPM
Oxygen Saturation	1-5 LPM 93 ± 3%
Oxygen Monitor	- Low : 82 % - Very Low : < 70 %
Weight	15.1 Kgs

