

# OXTM

# Venus 5

## PORTABLE OXYGEN CONCENTRATOR



**2.2 KG**



**5 Hours**



### Battery Duration

Setting	Mins
1	300
2	230
3	195
4	140
5	110

### Settings and bolus size (ml)

BPM	1	2	3	4	5
15	11	22	33	44	55
20	11	22	33	44	50
25	11	22	33	40	40
30	11	22	33	33	33
35	11	22	30	30	30
40	11	22	25	25	25

### Specification

<b>Oxygen Concentration</b>	87%-96% at all Flow Rates Based on atmospheric of 14.7 psoa (101 kPa) at 21 °C
<b>Setting</b>	1 ~ 5 Settings
<b>Power Requirements</b>	<u>AC power</u> 100~240 V AC (±10 %) 50/60 Hz in , 16.8 V DC, 6 A out <u>DC power</u> 12.0~16.8 V DC in, 16.8 V DC out
<b>Sound Level</b>	45 dBA at setting 2 (when measured at 1 m from front of device )
<b>Alarm Type</b>	◆ Low oxygen purity    ◆ Technical fault ◆ No breath            ◆ High breath rate ◆ Low battery           ◆ power failure ◆ No flow
<b>Outlet Pressure</b>	6~15 psi
<b>Dimension(mm)</b>	115(L) x 158 (W) x 255 (H) mm
<b>Weight</b>	4.84 LBS. 2.2 kgs (with internal battery)
<b>Battery Duration</b>	Up to 5 Hours
<b>Warranty</b>	3 years / 1 year for accessories
<b>Accessories</b>	Nasal cannula Power cord Carry bag Adaptor

## Oxytek Corp.

Email: [service@oxtm-o2.com](mailto:service@oxtm-o2.com) Tel: +886 4 2321 8700  
9F-2, No.239, Taiwan Blvd., West Dist., Taichung City 40350, Taiwan

